

Contact us

The Community Dental Service has a diaco dental chair located at the following clinic:

Middleton Community Clinic

Acre Close

Leeds

LS10 4HT

Tel: 0113 207 7064

If you would like to see the diaco dental chair in use please feel free to visit their website at: www.diaco.co.uk, or you can contact the clinic above for more information.

Help us get it right

The diaco dental chair is safe, simple and easy to use. It helps us ensure that our wheelchair and bariatric patients are as comfortable as can be and that they can receive the same treatment as all our other patients.

If you have a complaint, concern, comment or compliment, please let us know by speaking to a member of staff. We learn from your feedback and use the information you provide to improve and develop our services.

If you would like to talk to someone outside the service contact the **Patient Experience Team** on **0113 220 8585**, Monday to Friday 8.30am to 4.30pm or email **lch.pet@nhs.net**

We can make this information available in Braille, large print, audio or other languages on request.



Leeds Community
Healthcare
NHS Trust

Community Dental Service



Daco Dental Chair

Information on the Mobile
Dental Chair for wheelchair
and bariatric patients

The Diaco Dental Chair

The diaco dental chair is a fully mobile platform designed specifically for wheelchair and bariatric patients. The chair is suitable for manual and electric wheelchairs and will safely hold a maximum combined weight of 1100lbs/500kgs.



The benefits of the chair

The diaco chair allows us to increase the range of treatment we can offer to patients who are unable to use a standard dental chair. Our patients are able to remain in their own wheelchair at all times eliminating the need for manual handling, slings and hoisting. This not only makes it more comfortable but also safer and easier for patients and our clinical staff.



The fully adjustable head and back rest allows the patient to be reclined to the same position as a standard dental chair, this allows the dental team to work in the correct position ensuring that the patient is able to receive the best possible oral care.



Patient review

Following treatment in the diaco chair one patient commented:

“I backed up onto the platform in my wheelchair, applied my brakes and then had a headrest positioned behind my head. The platform which I was now on, still in my wheelchair, tilted backwards slowly and safely. I felt very secure. The best thing for me was being able to stay in my own wheelchair. I have to sit on a special gel cushion to ensure I have no pressure points on my skin, this is one thing that most paralysed people fear the most. Sitting on my own cushion in my own wheelchair is so reassuring, for me it eliminated this and so many other fears.

The tilting of the chair is smooth, it felt safe and I can see so many more disabled people now willing, and able, to visit the dentist.”